# RENASYS<sup>®</sup>-G STERILE

## Negative Pressure Wound Therapy Basic Gauze Dressing Kit O.R. Compatible

### Indications for Use

Smith & Nephew Wound Dressing Kits are intended to be used in conjunction with Smith & Nephew Negative Pressure Wound Therapy (NPWT) systems. NPWT systems are indicated for patients who would benefit from a suction device, particularly as the device may promote wound healing.

NPWT is appropriate for use on the following wounds:

- Pressure ulcers
- Diabetic/neuropathic ulcers
- Venous insufficiency ulcers
- Traumatic wounds
- Post-operative and dehisced surgical wounds
- Explored fistulas
- Skin flaps and grafts

### Contraindications

The use of the NPWT is contraindicated for:

- Untreated osteomyelitis
- Exposed arteries, veins or organs
- Eschar
- Malignancy in wound (with exception of palliative care to enhance quality of life)

### Warnings

1. Patients must be closely monitored for bleeding. If sudden or increased bleeding is observed, immediately discontinue therapy, take appropriate measures to stop bleeding and contact the treating clinician.
2. Patients suffering from difficult hemostasis or who are receiving anticoagulant therapy have an increased risk of bleeding. During therapy, avoid using hemostatic products that may increase the risk of bleeding, if disrupted.
3. Sharp edges or bone fragments in a wound must be covered or removed prior to using the NPWT system due to risk of puncturing organs or blood vessels while under negative pressure.
4. Do not use NPWT on exposed blood vessels or organs.
5. In the event defibrillation is required, disconnect the device from the wound dressing prior to defibrillation. Remove the wound dressing if it will interfere with defibrillation.
6. NPWT pumps are not MRI or CT compatible. Do not bring NPWT into the MRI suite or scanner range.
7. When operating, transporting, repairing or disposing of NPWT systems and accessories, the risk of infectious liquids being aspirated or contamination of the pump assembly through incorrect use cannot be eliminated. Universal precaution should be observed whenever working with potentially contaminated parts or equipment.
8. These devices have not been studied on pediatric patients. Patient size and weight should be considered when prescribing these devices.
9. NPWT is unsuitable for use in areas where there is danger of explosion (e.g., hyperbaric oxygen unit).
10. The canister kits are provided non-sterile and should not be used in a sterile field.

### Precautions

1. Precautions should be taken for patients who are or may be:
  - Receiving anticoagulant therapy or platelet aggregation inhibitors, actively bleeding or have weakened blood vessels or organs.
  - Suffering from difficult wound hemostasis.
  - Untreated for malnutrition.
  - Noncompliant or combative.
  - Suffering from wounds in close proximity to blood vessels or delicate fascia.

2. Infected wounds may require more frequent dressing changes. Regular monitoring of the wound should be maintained to check for signs of infection.
3. NPWT therapy should remain on in the CONSTANT mode for the duration of the treatment. There may be situations when the patient needs to be disconnected from the NPWT pump, for example for activities of daily living and diagnostic testing. If the patient needs to be disconnected the tubing should be clamped and the ends of the tubing protected. How long patients may be disconnected from the vacuum pump is a clinical decision based on individual characteristics of the patient and the wound. Factors to be considered would include; amount of drainage, location of the wound, integrity of the dressing seal, assessment of bacterial burden in the wound and patient's risk of infection.
4. As a condition of use, this device should only be used by qualified and authorized personnel. The user must have the necessary knowledge of the specific medical application for which NPWT is being used.
5. If the device has been at temperatures below freezing, the device must be brought to room temperature prior to use or the pumping unit may be damaged.
6. Ensure that tubing is installed completely and without any kinks to avoid leaks or blockages in the suction circuit.
7. Position the device and tubing appropriately to avoid the risks of causing a trip hazard and of the patient lying on the tubing.
8. When bathing/showering, the patient must clamp off the tubing and disconnect from the device.
9. Check the overflow protection/bacteria filter and replace as necessary or minimum of weekly with each canister change.
10. Underlying structures, such as tendons, ligaments and nerves should be covered with natural tissue or a non-adherent dressing layer.
11. The suction level should never be painful. If the patient reports discomfort with the suction level, consider reducing the pressure.
12. Always use the smallest canister volume possible – do not use the 800cc canister on patients with a high risk of bleeding.
13. Maintain regular monitoring of the device and wound site during therapy to ensure therapeutic treatment and patient comfort.
14. If any liquids penetrate the pump, discontinue use and return to your authorized provider for service.
15. As with all adhesive products apply and remove the dressing carefully from sensitive or fragile skin to avoid skin stripping, especially after frequent dressing changes.
16. Do not use if packaging is breached or damaged.

### Dressing Application

Use clean/aseptic or sterile techniques for application depending on institutional protocol. Only use the Smith & Nephew Dressing Kits approved for use with the system. If the objective is to bathe the wound surface for a period of time, the suction should be turned off during that period. When the irrigation port on the canister is not in use, the cap should be applied to the port.

#### Step 1


Cleanse wound bed and pat dry per protocol. Do not apply No-Sting Skin-Prep<sup>®</sup> directly onto open wound.


**Step 2** Apply skin sealant to periwound area.  
*Not included as part of Basic Sterile Kit.*

**Step 3** Trim single layer of non-adherent gauze and lay across wound bed.  
*Not included as part of Basic Sterile Kit.*

**Step 4** Cut the drain approximately 1" (2.5cm) shorter than the base of the wound.

**Step 5** Apply layer of saline-moistened antimicrobial gauze to wound bed. Position drain on top of gauze. For channel drain, wrap layer of gauze around drain.


*Caution: it is critical that the gauze is moistened with normal saline prior to filling the wound.*

**Step 6** Secure the drain in place.

**Step 7** Fill wound with saline-moistened antimicrobial gauze. Gauze should loosely fill to surface of wound. Avoid overpacking.  
*Saline not included as part of Basic Sterile Kit.*


**Step 8** Place transparent dressing over the wound and seal.

**Step 9** Secure the drain tubing with waterproof tape.

**Step 10** Connect the drain to the canister and begin suction, ensuring canister is installed correctly.

**Finished Dressing** will have a "raisin-like" appearance and be firm to the touch.  
*WARNING: The drain should never be placed in an unexplored fistula tract.*

### Dressing Changes

The Dressing Kit should be changed every 48 hours after the initial application of therapy. If no leak is present and the patient is comfortable, dressing changes should occur 2-3 times per week. In the event of heavy drainage or drainage with sediment, more frequent dressing changes may be needed.

Check dressings regularly. Infected wounds may require more frequent dressing changes. Monitoring of the wound must be maintained to check for signs of infection. Throughout treatment, monitor for any signs of local or systemic infection. If there are any signs of systemic infection or advancing infection at the wound site, contact the treating clinician immediately.

Made in USA for:

Smith & Nephew, Inc., St. Petersburg, FL 33716  
Smith & Nephew Medical Ltd.,  
Hull HU3 2BN, England

<sup>®</sup>Trademark of Smith & Nephew

© Smith & Nephew 2009 www.smith-nephew.com  
24/7 NPWT Clinical Support:  
1-866-998-NPWT (6798)

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**China**

Smith & Nephew Ltd.  
03-06,10 Floor,  
Headquarters Building,  
168 Xizang Road (M),  
HuangPu District,Shanghai,China  
200001  
Tel: +86-021-63503100  
Fax: +86-021-63501219  
info.china@smith-nephew.com

**Denmark**

Smith & Nephew A/S  
Slotsmarken 14  
2970 Hoersholm  
Danmark  
Tel: +45 4580 6100  
Fax: +45 4580 6151  
sn.dk@smith-nephew.com

**Finland**

Smith & Nephew Oy  
Äyritie 12 C  
01510 Vantaa  
Finland  
Tel: +358 (0)207 866 300  
Fax: +358 (0)207 866 320  
Tilaukset: +358 (0)207 866 333  
etunimi.sukunimi@smith-nephew.com

**France**

Smith & Nephew SAS  
Espace Novaxis  
25 Blvd. Marie et Alexandre Oyon  
72019 Le Mans Cedex 2  
France  
Tel: (33) 02 43 83 23 82  
(33) 02 43 83 24 14  
Fax: (33) 02 43 83 23 83  
scw@smith-nephew.com

**Germany**

Smith & Nephew GmbH  
Osterbrooksweg 71  
D-22869 Schenefeld  
Tel: +49 (040) 87 97 44 -0  
Fax: +49 (040) 87 97 44-375  
info@smith-nephew.com

**Hong Kong**

Smith & Nephew Ltd  
Unit 813-816, 8th Floor  
Delta House  
3 On Yiu Street  
Shatin NT  
Hong Kong  
Tel: +852 2648 7700  
Fax: +852 2648 7282  
hongkong@smith-nephew.com

**India**

Smith & Nephew Healthcare Pvt. Ltd.  
B-203, Dynasty Business Park,  
J. B. Nagar, Andheri Kurla Road,  
Andheri East,  
Mumbai - 400 059  
India  
Tel: +91 22 40055090  
Fax: +91 22 40055089  
india.admin@smith-nephew.com

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**Singapore**

Smith & Nephew Pte Ltd  
2 Jurong East Street 21  
IMM Building #04-32B  
Singapore 609601  
Tel: 65-62700552  
Fax: 65-62726698  
singapore.office2@smith-nephew.com

**South Africa**

Smith & Nephew (Pty) Ltd.  
PO Box 92  
Pinetown 3600  
South Africa  
Tel: +27 (0)31 242 8111  
Fax: +27 (0)31 242 8120  
lorraine.speed@smith-nephew.com

**South Korea**

Smith & Nephew Ltd.  
ASEM Tower(13F)  
Samsung-dong  
Gangnam-gu  
Seoul  
(Seoul 135-798 Korea)  
Tel: (02)6001-7575  
Fax: (02)6001-7560

**Spain**

Smith & Nephew S.A.  
Frutuós Gelabert, 2 y 4  
08970 Sant Joan Despí  
Barcelona  
España  
Tel: +34 93 373 7301  
Fax: +34 93 373 7453  
maria.moreno@smith-nephew.com

**Sweden**

Smith & Nephew AB  
Kråketorpsgatan 20  
Box 143  
431 22 Mölndal  
Sverige  
Tel: 031-746 58 00  
Fax: 031-87 05 32

**Switzerland**

Smith & Nephew AG  
Glutz Blotzheim-Strasse 1  
4502 Solothurn  
Switzerland  
Tel: +41 32 624 56 60  
Fax: +41 32 624 56 81  
swiss.info@smith-nephew.com

**Taiwan**

Smith & Nephew Taiwan  
9F - 3, No 130, Sec. 2  
Jhongsiao East Road  
Taipei City, Taiwan (R.O.C.)  
Tel: 886 2 3322 5020  
Fax: 886 2 3322 5021

**Thailand**

Smith & Nephew Ltd  
344/3 Soi Rongrien Yepun  
Rama IX Road  
Kwang Bangapi  
Khet Huay Kwang  
Bangkok 10320  
Thailand  
Tel: +66 2 719 6221-27  
Fax: +66 2 719 6220  
snlthai@loxinfo.co.th

**Turkey**

Baral Tibbi Malzemeler Ltd Sti  
Nilgun Sok. No. 12/12 06680  
Kavaklidere  
Ankara  
Turkey  
Tel: +90 312 427 1167  
bahadiralev@baral.com.tr

**United Arab Emirates**

Smith & Nephew FZE  
Level 4,  
National Bank of Abu Dhabi Bldg.  
Dubai Health Care City  
P.O.Box 9715  
Dubai  
United Arab Emirates  
Tel: +971 4 4299111  
Fax: +971 4 4299100  
snfze@emirates.net.ae

**United Kingdom**

Smith & Nephew Healthcare Ltd.  
Healthcare House  
Goulton Street  
Hull, HU3 4DJ  
United Kingdom  
Tel: 0800.590173  
Fax: 01482.222211  
advice@smith-nephew.com  
orders.uk@smith-nephew.com

**United States**

Smith & Nephew, Inc.  
Wound Management Division  
970 Lake Carillon Dr.  
St. Petersburg, FL 33716  
Tel: 1 800-876-1261  
Fax: 727-392-6914  
customercare.largo@smith-nephew.com  
1-866-998-NPWT (6798)