RENASYS[®]

High Output/Fistula Kit

Negative Pressure Wound Therapy

NOTE: The High Output Fistula Kit is not compatible with the RENASYS GO device.

Rx only

STERILE EO

2

Indications for Use

The RENASYS High Output Fistula Kit is intended to be used in conjunction with Smith & Nephew Negative Pressure Wound Therapy (NPWT) Systems to deliver negative pressure to the wound. Smith & Nephew NPWT Systems are indicated for patients who would benefit from a suction device, particularly as the device may promote wound healing.

NPWT is appropriate for use on the following wounds:

- Chronic
- Acute
- Traumatic
- Sub-acute and dehisced wounds
- Ulcers (such as pressure or diabetic)
- Partial-thickness burns
- Flaps and grafts

Contraindications

The use of NPWT is contraindicated for:

- Untreated osteomyelitis
- Exposed arteries, veins, organs or nerves
- Necrotic tissue with eschar present
- Malignancy in wound (with exception of palliative care to enhance quality of life)
- Non-enteric and unexplored fistulas
- Anastomotic sites

Warnings

1. Patients must be closely monitored for bleeding. If sudden or increased bleeding is observed, immediately discontinue therapy, take appropriate measures to stop bleeding and contact the treating clinician.
2. Patients suffering from difficult hemostasis or who are receiving anticoagulant therapy have an increased risk of bleeding. During therapy, avoid using hemostatic products that, if disrupted, may increase the risk of bleeding.
3. Sharp edges or bone fragments in a wound must be covered or

removed prior to using NPWT due to risk of puncturing organs or blood vessels while under negative pressure.

4. Do not tightly pack or force gauze into any areas of the wound.
5. Do not use NPWT on exposed blood vessels or organs.
6. In the event defibrillation is required, disconnect the NPWT device from the wound dressing prior to defibrillation. Remove the wound dressing if it will interfere with defibrillation.
7. The NPWT device is not MRI or CT compatible. Do not bring into the MRI suite or scanner range.

8. When operating, transporting, repairing or disposing of NPWT devices and accessories, the risk of infectious liquids being aspirated, or contamination of the device assembly through incorrect use, cannot be eliminated. Universal precaution should be observed whenever working with potentially contaminated parts or equipment.

9. NPWT has not been studied on pediatric patients. Patient size and weight should be considered when prescribing NPWT devices.
10. NPWT is unsuitable for use in areas where there is danger of explosion (e.g., hyperbaric oxygen unit).

Precautions

1. Precautions should be taken for patients who are or may be:
 - Receiving anticoagulant therapy or platelet aggregation inhibitors, actively bleeding or have weakened blood vessels or organs.
 - Suffering from difficult wound hemostasis.
 - Untreated for malnutrition.
 - Noncompliant or combative.
 - Suffering from wounds in close proximity to blood vessels or delicate fascia.
2. The use of NPWT presents a risk of tissue in-growth. Tissue in-growth may be reduced by reducing therapy pressure, using a wound contact layer or by increasing the frequency of dressing changes.
3. Infected wounds may require more frequent dressing changes. Regular monitoring of the wound should be maintained to check for signs of infection.
4. NPWT should remain on in the CONTINUOUS mode for the duration of the treatment. There may be situations when the patient needs to be disconnected from the NPWT device, for example for activities of daily living and diagnostic testing. If the patient

needs to be disconnected, the dressing tubing should be clamped and the ends of the tubing protected. How long patients may be disconnected from the NPWT device is a clinical decision based on individual characteristics of the patient and the wound. Factors to be considered would include; amount of drainage, location of the wound, integrity of the dressing seal, assessment of bacterial burden in the wound and patient's risk of infection.

5. As a condition of use, the NPWT device should only be used by qualified and authorized personnel. The user must have the necessary knowledge of the specific medical application for which NPWT is being used.
6. If the NPWT device has been at temperatures below freezing, the device must be brought to room temperature prior to use or the pumping unit may be damaged.
7. Ensure that tubing is installed completely and without any kinks to avoid leaks or blockages in the vacuum circuit.
8. Position the NPWT device and tubing appropriately to avoid the risks of causing a trip hazard and of the patient lying on the tubing.
9. When bathing/showering, the patient must clamp off the dressing tubing and disconnect from the NPWT device.

10. Inspect the overflow protection/bacteria filter on the canister and replace the canister as necessary. At minimum, the canister should be changed weekly.

11. Underlying structures, such as tendons, ligaments and nerves should be covered with natural tissue or a non-adherent dressing layer prior to applying the NPWT dressing kit.
12. NPWT should not be painful. If the patient reports discomfort, consider reducing the pressure.
13. Always use the smallest canister volume possible – do not use the 800ml canister on patients with a high risk of bleeding.
14. Maintain regular monitoring of the NPWT device and wound site during therapy to ensure therapeutic treatment and patient comfort.
15. If any liquids penetrate the NPWT device, discontinue use and return to your authorized provider for service.
16. As with all adhesive products apply and remove the dressing carefully from sensitive or fragile skin to avoid skin stripping, especially after frequent dressing changes.
17. Do not use if packaging is breached or damaged.
18. Canister kits are provided non-sterile and should not be used in a sterile field.

Precautions Specific to the RENASYS High Output Fistula Kit

Pressure Settings Recommendations

- Utilize device pressure settings between 40-120mmHg.
- The use of pressure settings which are higher than the recommended maximum of 120mmHg is a decision to be made by the treating physician/clinician and the individual needs of the patient should be taken into account.

Fistula Precautions

- Never lay drain directly in contact with the intestinal mucosa at the top of fistula opening.
- During the course of treatment patient's fluid levels must be closely monitored.

Recommendations Based on Fistula Location and Fluid Description

The anatomical location of the fistula in the intestine will influence the amount and type of the output:

- The higher the fistula occurs in the bowel, the thinner, more caustic and heavier the output.
- The lower the fistula occurs in the bowel, the thicker or pastier the output.

Dressing Application

Use aseptic or sterile techniques for application depending on institutional protocol. Use only Smith & Nephew Dressing Kits approved for use with the NPWT device.

NOTE: The drain in the High Output Fistula Kit is equipped with an irrigation port for irrigating the wound surface with an appropriate solution. If the objective is to bathe the wound surface for a period of time, the NPWT device should be turned off during that period. When the irrigation port is not in use, the cap should be applied to the port.

Step 1: Cleanse wound area and pat dry per protocol. Thorough wound cleansing should occur with each dressing change.


Step 2: Apply skin sealant to periwound area.
• Do not apply No-Sting Skin-Prep[®] directly on open wound.

Step 3: Cover wound bed, including fistula opening, with non-adherent gauze or with one layer of saline moistened gauze.

CAUTION: It is important that the gauze is moistened with normal saline prior to filling the wound.

High Output/Thin Consistency Dressing Procedure

Step 4: Place 28 Fr drain in close proximity to fistula opening and in dependent/inferior position.

Caution: Utilize stoma paste underneath drain at wound edge to minimize pressure point.

WARNING: The drain should never be placed directly in the fistula tract.

Step 5: Fill remaining wound deficit with additional saline moistened gauze.

Sediment Laden Dressing Procedure

Step 4: Place 28 Fr drain in close proximity to fistula opening and in dependent/inferior position. Fistula opening should remain exposed.

WARNING: The drain should never be placed directly in the fistula tract.

Step 5: Fill wound deficit with additional saline moistened gauze and continue to leave fistula opening exposed.

Modify the application of the High Output Fistula Kit as described below in accordance with the location and the fluid properties.

Step 6: Place transparent dressing over the wound and seal.

Step 7: Secure the 28 Fr drain with waterproof tape.

Step 8: Remove the canister tube supplied with the canister and discard. Attach the canister to the RENASYS device.

Step 9: Connect the High Output Fistula Drain directly to the canister by joining the blue connector of the drain to the canister input.

PRECAUTION: Ensure dressing tubing clamp is open and there are no kinks in the tubing.

Step 10: Activate the RENASYS device on CONTINUOUS mode beginning with 80mmHG. Check that the dressing has a good seal. Finished dressings should be firm to the touch. If required, adjust pressure setting to the desired level.

Dressing Changes

PRECAUTION: Before removing the dressing, clamp the dressing tubing immediately before turning off the device.

Step 1: Remove and dispose of dressing per institutional protocol. Thoroughly inspect the wound to ensure all gauze pieces have been removed.

Step 2: As with all adhesive products, apply and remove the dressing carefully from sensitive or fragile skin to avoid skin stripping, especially if frequent dressing changes will be needed.

Step 3: Dressings should be changed every 48 hours, or more frequently, based on continual monitoring of patient condition. In the event of heavy drainage, drainage with sediment or infected

wounds, more frequent dressing changes may be needed.

Step 4: Check dressings regularly and monitor the wound to check for signs of infection or bleeding. Infected wounds may require more frequent dressing changes. If there are any signs of systemic infection or advancing infection at the wound site, contact the treating clinician immediately.

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China

Smith & Nephew Ltd.
03-06,10 Floor,
Headquarters Building,
168 Xizang Road (M),
HuangPu District,Shanghai,China
200001
Tel: +86-021-63503100
Fax: +86-021-63501219
info.china@smith-nephew.com

Denmark

Smith & Nephew A/S
Slotsmarken 14
2970 Hoersholm
Danmark
Tel: +45 4580 6100
Fax: +45 4580 6151
sn.dk@smith-nephew.com

Finland

Smith & Nephew Oy
Äyritie 12 C
01510 Vantaa
Finland
Tel: +358 (0)207 866 300
Fax: +358 (0)207 866 320
Tilaukset: +358 (0)207 866 333
etunimi.sukunimi@smith-nephew.com

France

Smith & Nephew SAS
Espace Novaxis
25 Blvd. Marie et Alexandre Oyon
72019 Le Mans Cedex 2
France
Tel: (33) 02 43 83 23 82
(33) 02 43 83 24 14
Fax: (33) 02 43 83 23 83
scw@smith-nephew.com

Germany

Smith & Nephew GmbH
Osterbrooksweg 71
D-22869 Schenefeld
Tel: +49 (040) 87 97 44 -0
Fax: +49 (040) 87 97 44-375
info@smith-nephew.com

Hong Kong

Smith & Nephew Ltd
Unit 813-816, 8th Floor
Delta House
3 On Yiu Street
Shatin NT
Hong Kong
Tel: +852 2648 7700
Fax: +852 2648 7282
hongkong@smith-nephew.com

India

Smith & Nephew Healthcare Pvt. Ltd.
B-203, Dynasty Business Park,
J. B. Nagar, Andheri Kurla Road,
Andheri East,
Mumbai - 400 059
India
Tel: +91 22 40055090
Fax: +91 22 40055089
india.admin@smith-nephew.com

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Singapore

Smith & Nephew Pte Ltd
2 Jurong East Street 21
IMM Building #04-32B
Singapore 609601
Tel: 65-62700552
Fax: 65-62726698
singapore.office2@smith-nephew.com

South Africa

Smith & Nephew (Pty) Ltd.
PO Box 92
Pinetown 3600
South Africa
Tel: +27 (0)31 242 8111
Fax: +27 (0)31 242 8120
lorraine.speed@smith-nephew.com

South Korea

Smith & Nephew Ltd.
ASEM Tower(13F)
Samsung-dong
Gangnam-gu
Seoul
(Seoul 135-798 Korea)
Tel: (02)6001-7575
Fax: (02)6001-7560

Spain

Smith & Nephew S.A.
Frutuós Gelabert, 2 y 4
08970 Sant Joan Despí
Barcelona
España
Tel: +34 93 373 7301
Fax: +34 93 373 7453
maria.moreno@smith-nephew.com

Sweden

Smith & Nephew AB
Kråketorpsgatan 20
Box 143
431 22 Mölndal
Sverige
Tel: 031-746 58 00
Fax: 031-87 05 32

Switzerland

Smith & Nephew AG
Glutz Blotzheim-Strasse 1
4502 Solothurn
Switzerland
Tel: +41 32 624 56 60
Fax: +41 32 624 56 81
swiss.info@smith-nephew.com

Taiwan

Smith & Nephew Taiwan
9F-3, No. 130, Sec. 2
Jhongsiao East Road
Taipei City, Taiwan (R.O.C.)
Tel: 886 2 3322 5020
Fax: 886 2 3322 5021

Thailand

Smith & Nephew Ltd
344/3 Soi Rongrien Yepun
Rama IX Road
Kwang Bangapi
Khet Huay Kwang
Bangkok 10320
Thailand
Tel: +66 2 719 6221-27
Fax: +66 2 719 6220
snlthai@loxinfo.co.th

Turkey

Baral Tibbi Malzemeler Ltd Sti
Nilgun Sok. No.12/12 06680
Kavaklıdere
Ankara
Turkey
Tel: +90 312 427 1167
bahadiralev@baral.com.tr

United Arab Emirates

Smith & Nephew FZE
Level 4,
National Bank of Abu Dhabi Bldg.
Dubai Health Care City
P.O.Box 9715
Dubai
United Arab Emirates
Tel: +971 4 4299111
Fax: +971 4 4299100
snfze@emirates.net.ae

United Kingdom

Smith & Nephew Healthcare Ltd.
Healthcare House
Goulton Street
Hull, HU3 4DJ
United Kingdom
Tel: 0800.590173
Fax: 01482.222211
advice@smith-nephew.com
orders.uk@smith-nephew.com

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 Smith & Nephew Medical Ltd.,
Hull HU3 2BN, England

®Trademark of Smith & Nephew
© Smith & Nephew 2011 www.smith-nephew.com

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Smith & Nephew, Inc.
Smith & Nephew, Inc., St. Petersburg, FL 33716