



RENASYS-F/AB  
Abdominal Dressing Kit

 **smith&nephew**  
**RENASYS<sup>®</sup>-F/AB**  
Abdominal Dressing Kit

# Your solution for managing open abdominal wounds

Patients with critical open abdominal wounds require specialized treatment. When primary closure is not possible or repeat entries are necessary, you may require a comprehensive Negative Pressure Wound Therapy (NPWT) solution designed specifically for the open abdomen.

The RENASYS<sup>®</sup>-F/AB Abdominal Dressing Kit from Smith & Nephew is an intuitive NPWT solution for temporary abdominal closure (TAC). This abdominal kit can be used with the RENASYS EZ Plus device, a powerful NPWT pump that is compatible with all RENASYS Dressing Kits. This improves efficiencies by eliminating the need for a dedicated abdominal device. In addition, you have the flexibility to leverage your RENASYS equipment across multiple care settings, patients and wound types.

## Why choose the RENASYS-F/AB?

### A cost-effective solution


The RENASYS-F/AB Dressing Kit is recommended for use with the powerful and intuitive RENASYS EZ Plus device, which can also be used with the RENASYS-F/P Foam and RENASYS-G/P Gauze Dressing Kits. This simplifies supply management and reduces training.

### Fast, simple application

This RENASYS-F/AB Kit includes everything required to achieve temporary abdominal closure. The organ protection layer is quick and easy to apply, saving costly OR time. Two pieces of pre-shaped foam allow you to meet specific needs and circumstances of the patient's open abdomen management. And the port dressing is simple and familiar.

### Supports clinical goals

The RENASYS-F/AB Kit protects and controls the abdominal contents by providing separation between the abdominal wall and the viscera. It also removes exudate and allows estimation of fluid loss, helps reduce edema, provides medial tension to slow fascial retraction, and supports timely primary fascial closure of the open abdomen.



## Flexible and complete

The RENASYS-F/AB Kit is part of the complete RENASYS solution for NPWT. With a portfolio of compatible dressing kits offering a choice of foam and gauze interfaces, devices and accessories, you can cost-effectively treat a wide variety of chronic, acute and abdominal wounds.

For critical abdominal wounds, the RENASYS-F/AB Kit works effectively with the RENASYS EZ Plus device, which offers simple controls, a high flow rate, and multiple pressure settings to help you customize the delivery of NPWT for high volume output wounds.



## RENASYS<sup>◇</sup>-F/AB: Ordering information

### Dressing Kits

Product Description	Product Code
<b>RENASYS-F/AB</b> Dressing Kit includes: Organ Protection Layer (OPL) x 1 Perforated Pre-shaped Foam x 2 Transparent Film x 6 Port Dressing x 1	<b>Smith &amp; Nephew</b> 66800713  <b>Acute Care (UHS)</b> 2010175



### Compatible Devices

<b>RENASYS EZ Plus</b>	<b>Smith &amp; Nephew</b> 66800697  <b>Acute Care (UHS)</b> ZNP
------------------------	---



### Compatible Canisters

<b>RENASYS 250ml Canister</b>	<b>Smith &amp; Nephew</b> 66800058  <b>Acute Care (UHS)</b> 2010150
-------------------------------	---



<b>RENASYS 800ml Canister</b>	<b>Smith &amp; Nephew</b> 66800423  <b>Acute Care (UHS)</b> 2009250
-------------------------------	---



## Clinical and support services

Smith & Nephew is committed to supporting caregivers and their patients through the entire continuum of care. We partner with national acute and home care distributors to provide a wide range of clinical services, so you have access to the expertise and resources you and your patients need at every point of care.

### Product and clinical support

- Product instructions and clinical application techniques available on website and DVD
- Product manuals, NPWT clinical guidelines and quick reference guides
- Dedicated team of home care specialists
- Clinical in-servicing and continued product training
- Patient education and support

### Customer and patient hotlines

- 24/7 clinical and patient hotlines
- Customer Care and Technical Support call centers

### Logistics and fleet management

- Patient equipment delivery, set up, service and maintenance for acute and home care
- Third-party billing management

In addition, we offer education and information about the therapeutic area of wound management. From interactive presentations and education modules to peer learning opportunities, Smith & Nephew provides ongoing professional development tools to expand your understanding of wound care trends and treatments.

### Clinical education

- Peer-to-peer learning through our National Science Summit Program
- Accredited online education modules
- Educational videos and webcasts

With over 150 years of experience in advanced wound care, Smith & Nephew is an industry leader providing innovative solutions to meet the needs of chronic, acute and traumatic wounds across all care settings.

#### Wound Management

Smith & Nephew, Inc.

**USA**  
970 Lake Carillon Drive  
Suite 110  
St. Petersburg, FL 33716

Customer Care Center  
1-800-876-1261  
T 727-392-1261  
F 727-392-6914

**Canada**  
2250 Alfred-Nobel Blvd.  
Suite 300  
St-Laurent, Quebec  
H4S 2C9

Customer Action Center  
T 1-800-463-7439  
F 1-800-671-9140

[www.smith-nephew.com](http://www.smith-nephew.com)  
[www.myrenasys.com](http://www.myrenasys.com)